

Product information

VistaScan image plate cleaning wipe



VistaScan at a glance

- Image plate cleaning wipe with excellent material compatibility.
- Specially formulated for the VistaScan image plate.
- Protects the value of dental equipment thanks to its tested material compatibility.

Properties

The VistaScan Image Plate Cleaning Wipe is a special intensive cleaning wipe, particularly suitable for the VistaScan Image Plate.

Product composition

VistaScan image plate cleaning wipes contain 4 g 2-propanol.

Directions for use

Use only the VistaScan image plate cleaning wipe for VistaScan image plates. Carefully wipe off the image plate with the cleaning wipe and allow to dry. Only clean undamaged image plates.

Environmental impact

The active ingredient solution of the VistaScan image plate cleaning wipe is readily biodegradable. The packing is made of polyethylene and can therefore be exploited both materially and thermally.

Shelf life

4 years

Container size

10 wipes per folding box

Storage

Store preparation in a cool place, but not below 5°C.

Hazard warnings

VistaScan image plate cleaning wipes are classified and labelled according to the CLP Regulation: see product label and safety data sheet.

Independent expert opinions – in-house investigations

The reports are available on request.



The recommendations concerning our products are based on exhaustive in-house and external scientific investigations and are provided to the best of our knowledge. Our liability is limited here based on our „General Sales and Delivery Terms“ that are available upon request at any time.

Sales

Dürr Dental SE
Höpfigheimer Str. 17
74321 Bietigheim-Bissingen
Germany
www.duerrdental.com
info@duerrdental.com

Manufacturer

orochemie GmbH + Co. KG
Max-Planck-Straße 27
70806 Kornwestheim
Germany
www.orochemie.de
info@orochemie.de

